



SEM


SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com







Tomosynthesis Phantom NS

Item No. PH-56
Weight 2kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Allowing evaluation of reconstruction slices and uniformity in the measurement of slice thickness through showing the images numerically and graphically

Training skills / Applications

Set includes

1.Verification of reconstruction interval / 2.Slice thickness / 3.Uniformity
1 reconstruction positioning unit / 1 slice thickness unit / 1 uniformity unit / 1 height setting rack / manual

Size (approx.)
Packing size (approx.)
Packing weight (approx.)
Materials
APPLICATIONS
Update

W7 x D15 x H25 cm W2.7 x D6 x H9.8 in.
W46 x D31 x H17cm / W18.1 x D12.2 x H6.7in
2 kg / 4.4 lb
Acrylic resin, bakelite, aluminum, copper, stainless
Tomosynthesis
August 10, 2020