



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




Tough Bone Phantom BE-T, BE-H, BE-NWD

Item No. PH-41
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

Human bone substitute phantoms to simulate body structure in combination with PH-40 and PH-42

Features

Size (approx.)

Materials

A stable, high quality and shatter-free phantom for radiotherapy.

Refer to catalog pdf

Epoxy resin

VARIATIONS

Dosimetry cavities

Update

BE-T-2005 Compact Bone 200x200x5mm/ 8x8x0.2 in / BE-T-2010 Compact Bone 200x200x10mm/ 8x8x0.4 in / BE-T-2020 Compact Bone 200x200x20mm/ 8x8x 0.8 in / BE-H-2005 Cortical Bone 200x200x5 mm/ 8x8x0.2 in / BE-H-2010 Cortical Bone 200x200x10mm/ 8x8x0.4 in / BE-H-2020 Cortical Bone 200x200x20mm / 8x8x0.8 in / BE-N-2005 Inner Bone 200x200x5 mm/ 8x8x0.2 in / BE-N-2010 Inner Bone 200x200x10 mm/ 8x8x0.4 in / BE-N-2020 Inner Bone 200x 200x 20 mm/ 8x 8x 0.8 in / BE-H-3005 Cortical Bone 300x 300x 5 mm/ 12x 12x 0.2 in / BE-H-3010 Cortical Bone 300x 300x 10mm/12x 12x 0.4 in / BE-H-3020 Cortical Bone 300x 300x 20 mm/ 12x 12x 0.8 in / BE-N-3005 Inner Bone 300 x 300 x 5 mm / 12 x 12 x 0.2 in / BE-N-3010 Inner Bone 300 x 300 x 10 mm / 12 x 12 x 0.4 in / BE-N-3020 Inner Bone 300 x 300 x 20 mm / 12 x 12 x 0.8 in
Tough series phantoms can be ordered with cavities and plugs. Specify your chamber's manufacturer and model number. Let us have dimensional drawings of the chambers you are using to estimate cost.
July 3, 2020