



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com




MRI Quality Assurance Phantom MHR

Item No. PH-31
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Radiology Phantom



Product Description

This QA phantom for MRI allows to evaluate the slice thickness, spatial resolution, uniformity and geometric distortion as well as contrast

Features

1. Uniformity is maintained under the high magnetic field of 3.0 Tesla
2. Uniformity provides high precision evaluation for other parameters

Set includes

Size (approx.)

Packing size (approx.)

Materials

Update

1 phantom unit A / 1 phantom unit B / 1 liquid paraffin / 1 dropper/ 5 NiCl 50ml (5, 10, 15, 20, 25 mmol) / 7 sample bottle (13.5ml) / 1 funnel / 1 petroleum jelly / 1 screwdriver / 1 extra screws / 1 storage case / manual
22 dia. x 14(H)cm x 2 types 8.7 dia. x 5.5 in
W66xD54xH34 cm / W26xD21.3xH13.4 inch
acrylic resin MRI contrast solution: nickel dichloride (NiCl)
July 9, 2020