



SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Ultrasound Neonatal Head Phantom (Abnormal type)

Item No. US-14b
Weight 5 kg
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Ultrasound Phantom



Product Description

World's first ultrasound neonatal head phantom with hydrocephalus

Features

1. Accurate depiction of a newborn's cerebral anatomy
2. Realistic user experience with its life-like soft touch.

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Related products

Production ? Development Supervision

ANATOMY

Update

Scanning of brain anatomy in sagittal (angled parasagittal), coronal and transverse planes via any fontanel / Palpation of anterior, posterior and mastoid fontanel
1 ultrasound phantom / 1 pillow / 1 storage case / 1 talcum powder/ 1 instruction manual
W13 x D14 x H20cm W5.1 x D5.5 x H7.9in
W35 x D25 x H19cm / W13.8 x D9.8 x H7.5in
1.5kg / 3.3lb
5kg / 11lbs
Polyurethane elastomer, Latex free
US-14a Ultrasound Neonatal Head Phantom (Normal type) US-13 Infant Hip Sonography Training Phantom Newborn Whole Body Phantom "PBU-80"
Charlotte Henningsen ,MS,RT(R),RDMSRVT,FSDMS,FAIUM Chair & Professor - Sonography Department Adventist University of Health Sciences
skull, anterior fontanel, posterior fontanel, hypertrophied lateral ventricle
August 9, 2020