



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com

Ultrasound-Guided Thoracentesis/ Pericardiocentesis Simulator

Item No. MW17
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories:Advanced Trauma Life Support (ATLS)



Product Description

Simulate thoracentesis on both a torso and SP / Practice patient safety during pericardiocentesis with ultrasound guidance

Features

1. Palpable ribs 2. Realistic needle-tip insertion and resistance 3. Both thoracentesis and pericardiocentesis in one unit

Training skills / Applications

Set includes

Size (approx.)

Packing size (approx.)

Weight (approx.)

Packing weight (approx.)

Materials

Replacement parts

Related products

Production ? Development Supervision

Recommended Devices

Notes

Update

Patient positioning / Recognition of anatomical landmarks by ultrasound / Visualization of pericardial fluid under ultrasound scanning / Landmark palpation / Assessment of level and volume of pleural effusion / Determination of insertion site / Needle insertion and collection of fluid

1 adult chest model / 1 mid-axially line unit / 1 mid-scapular line unit / 1 pericardiocentesis unit / 1 pillow for positioning / 1 explanation model for thoracentesis / 1 instruction manual / 1 irrigator / 1 funnel / 2 syringes (50ml) / 1 joint hose / 1 tube with three-way stopcock / 1 plastic jar

Torso size: W38 x D25 x H48 cm , 15 x 9.9 x 18.9 in /
Thoracentesis pad size: W16 x D7 x H21 cm , 6.3 x 2.76 x 8.3 in / Perocardio pad size: W16 x D14 x H21 cm , 6.3 x 5.6 x 8.3 in

W66 x D49 x H46 cm / W26 x D19 x H18 in

9.65kg / 4.4lb

12kg / 26.5lb

Soft resin / Latex free

11383-010 mid axially line unit puncture pad / 11383-020 mid-scapular line unit puncture pads / 11383-030 replacement lung / 11394-010 pad for MW15

[MW4 Ultrasound-Guided Thoracentesis Simulator M4A](#)
[Ultrasound-Guided Thoracentesis Simulator - Strap-on set - MW15 Ultrasound-Guided Pericardiocentesis Simulator](#)

Takahiro Amano, M.D., Vice President, Senior Vice-Dean, Postgraduate School, Professor and Director, Center of Postgraduate Medical Education, International University of Health and Welfare Honorary Director, Sanno Medical Center

Thoracentesis: 22 G - 23 G needle / Pericardiocentesis: 18 G needle

Needles are not included.

October 20, 2020