



**SEM**

## SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

**Mobile :** +91 88495 63724  
**Mobile 1:** +91 98791 03905

**Email:** sem@semtrainers.com  
**Website:** www.semtrainers.com

## Ultrasound-Guided Thoracentesis Simulator

Item No. MW4  
Weight  
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

**SKU:**

**Categories:**Ultrasound Phantom



### Product Description

Simulate thoracentesis on both a torso and SP / Ultrasound-Guided Thoracentesis Simulator features two types of puncture units: mid-scapular line access and mid-axillary line access. Strap-on (wearable) puncture units facilitates hybrid training sessions with simulated patients.

## Features

## Training skills / Applications

## Set includes

## Size (approx.)

## Packing size (approx.)

## Weight (approx.)

## Packing weight (approx.)

## Materials

## Replacement parts

## Optional parts

## Related products

## Production ? Development Supervision

## Recommended Devices

## Notes

## Update

1. Strap-on units to practice patient positioning and communication 2. Two sites for access: right mid-scapular line and left mid-axillary line 3. Volume of pleural effusion can be controlled to set different levels of difficulty.

Patient positioning / Recognition of anatomical landmarks by ultrasound / Assessment of level and volume of pleural effusion / Determination of insertion site / Needle insertion and collection of fluid

1 adult chest model / 1 mid-axially line unit / 1 mid-scapular line unit / 1 pillow for positioning / 1 explanation model for thoracentesis / 1 instruction manual / 1 irrigator / 1 funnel / 1 syringes (50ml) / 1 tube with three-way stopcock / 1 plastic jar

Torso size: W40 x D25 x H50 cm , 15.7 x 9.9 x 19.7 in / Thoracentesis pad size: W16 x D7 x H21 cm , 6.3 x 2.76 x 8.3 in

W65 x D48 x H48 cm / W25.5 x D19 x H19 in

9.65kg / 4.4lb

12kg / 26.5lb

Soft resin / Latex free

11383-010 mid axially line unit puncture pad / 11383-020 mid-scapular line unit puncture pads / 11383-030 replacement lung

11394-010 pad for MW15

[MW17 Ultrasound-Guided Thoracentesis/ Pericardiocentesis Simulator](#) [M4A Ultrasound-Guided Thoracentesis Simulator - Strap-on set](#) - [MW15 Ultrasound-Guided Pericardiocentesis Simulator](#)

Takahiro Amano, M.D., Vice President, Senior Vice-Dean, Postgraduate School, Professor and Director, Center of Postgraduate Medical Education, International University of Health and Welfare Honorary Director, Sanno Medical Center

22 G - 23 G needle

Needles are not included.

October 20, 2020