



SEM

SEM TRAINERS & SYSTEMS

WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE

SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile : +91 88495 63724
Mobile 1: +91 98791 03905

Email: sem@semtrainers.com
Website: www.semtrainers.com


Lumbar Puncture Simulator IIA

Item No. M43B
Weight
Brand KYOTO KAGAKU, JAPAN

[Read More](#)

SKU:

Categories: Epidural and Lumbar



Product Description

World's best training simulator for lumbar puncture and epidural anesthesia / New stable supporting stand for one person and team training scenarios/ Five types of puncture blocks enhance training with different types of challenges

Features

Training skills / Applications

Set includes

Size (approx.)

Materials

Replacement parts

Optional parts

Related products

Production ? Development Supervision

Recommended Devices

Notes

Update

1. Realistic skin and tissue resistance during insertion
Measure CSF fluid pressure, collect CSF fluid and practice epidural anesthesia procedures
2. Simulates administration of epidural nerve block with realistic needle-tip sensation
3. Simulates loss of resistance technique

Palpation of landmarks / Lumbar puncture / CSF collection / CSF pressure measurement / Epidural anesthesia / Anatomical understanding

1 lumbar region model / 6 lumbar puncture blocks: 2 normal CSF, 1 obesity CSF, 1 senior CSF, 1 senior obesity CSF, 1 epidural / 1 lumbar region skin cover / 1 supporting stand / 1 support base for team teaching / 1 lumbar spine model / 1 irrigation bag / 1 syringe / 1 stand for irrigation bag / 1 instruction manual / 1 storage case

W33 x D21 x H30 cm / W13 x D8.3 x H11.8 in

Soft resin, hard resin / Latex free

11348-090 1 LP block (normal CSF)/ 11348-110 1 LP block (obesity CSF)/ 11348-120 1 LP block (senior CSF)/ 11348-130 1 LP block (senior obesity CSF)/ 11348-140 1 epidural block/ 11348-150 1 lumbar region skin cover light tone/11348-150-D 1 lumbar region skin cover dark tone /11348-160 supporting stand (two ways)

11348-190 ultrasound lumbar puncture-epidural block / 11348-230 ultrasound lumbar region skin cover for M43E

[MW3 Epidural Anesthesia Simulator](#) [M43D Pediatric Lumbar Puncture Simulator II](#) [M43E Ultrasound Compatible Lumbar Puncture/ Epidural Simulator](#)

Takahiro Amano, M.D., Vice President, Senior Vice-Dean, Postgraduate School, Professor and Director, Center of Postgraduate Medical Education, International University of Health and Welfare Honorary Director, Sanno Medical Center

Epidural needle: 18G / Spinal needle: 21G

Needles are not included.

October 1, 2020