
SEM TRAINERS & SYSTEMS
WHERE TIME MOVES AHEAD TO KEEP PACE WITH KNOWLEDGE
SEM- Scientific Educative Methods in Science, Engineering & Medicine

Mobile :	+91 88495 63724	**Email:**	sem@semtrainers.com
Mobile 1:	+91 98791 03905	**Website:**	www.semtrainers.com



Outer chest skin for Deluxe CRiSis™ trunk simulator

Item No.	1019763
MPN:	LF03191U

Read More

SKU:

Categories:ALS Adult

Product Description

Outer chest skin for Deluxe CRiSis™ trunk simulator W44135

Why Choose the Life/form® Outer Chest Skin for Deluxe CRiSis™ Trunk Simulator for Your Simulation Lab in India?

The Life/form® Outer Chest Skin for Deluxe CRiSis™ Trunk Simulator (Item No. 1019763, MPN: LF03191U) by Nasco Healthcare, USA is the genuine OEM replacement outer chest skin for the Deluxe CRiSis™ trunk simulator. The outer chest skin is the anatomically detailed surface layer that trainees interact with directly — it is what their hands contact during CPR hand placement and compression, where defibrillation pads are positioned, and what defines the visual and tactile realism of the chest for assessment training. It is a distinct component from the Trunk Shell (LF03190U), which is

the underlying structural body of the torso.

The outer chest skin degrades differently from the trunk shell — it is subject to direct mechanical wear from repeated CPR compressions, defibrillation pad contact, and general handling across training sessions. As it degrades, the surface texture changes, the anatomical landmarks become less defined, and the material may crack or lose elasticity — all of which reduce the training fidelity that trainees rely on for correct hand placement, landmark identification, and compression technique development. Replacing the outer chest skin with the genuine OEM LF03191U part restores the correct surface texture, correct anatomical detail, and correct material compliance that CPR and defibrillation training require. Only the genuine OEM skin is guaranteed to fit the Deluxe CRiSis™ trunk frame with the correct dimensions, correct mounting geometry, and correct material specification. For high-throughput simulation labs running BLS, ACLS, and ALS certification programmes on the CRiSis™ platform, the outer chest skin is a consumable with a finite service life that should be stocked as a scheduled replacement item.

Why Buy from SEM Trainers? SEM Trainers & Systems is an authorized distributor of Nasco Healthcare products in India with over 30 years of experience supplying medical simulators and genuine OEM spare parts to 1,500+ institutions across 28 states. We provide pan-India delivery and dedicated after-sales support. All products carry genuine manufacturer warranty.

Frequently Asked Questions

Q: What is the difference between the Outer Chest Skin (LF03191U) and the Trunk Shell (LF03190U)? A: The outer chest skin (LF03191U) is the external anatomical surface layer — the soft skin that trainees interact with during CPR and defibrillation training. The trunk shell (LF03190U) is the underlying rigid structural body of the torso. Both are separate consumable/replacement parts for the Deluxe CRiSis™ trunk simulator.

Q: Which CRiSis™ simulators is this outer chest skin compatible with? A: LF03191U is designed for the Deluxe CRiSis™ trunk simulator by Nasco Healthcare. Contact us to confirm compatibility with your specific configuration before ordering.

Q: How do I order and get support from SEM Trainers? A: Contact us at sem@semtrainers.com or call +91-88495 63724. We supply genuine Nasco Healthcare spare parts with pan-India delivery.